



**Headquarters**

#951, Woram-Dong, Dalseo-Gu, Daegu, Korea / Tel. +82-53-583-2804~5 Fax. +82-53-583-2806 [www.dentis.co.kr](http://www.dentis.co.kr)

**Sales Head Office(Seoul)**

#907-908, Leaders Tower 60 -15, Gasan-Dong Geumcheon-Gu Seoul, Korea

**DENTIS USA**

5455 Wilshire Blvd #950 LA CA 90036 / Tel. +1-323-677-4363~5 Fax. +1-323-677-4366

**DENTIS TAIWAN**

11F-2, No.65, Songde Rd., Xinyi Dist., Taipei City 110, Taiwan / Tel : + 886-2-2729-6200

Issued date : 1<sup>st</sup> of December, 2011 ©2011 DENTIS CO.,LTD.

Version : Rev.201112-1A All Right reserved.

Specifications are subject to change without notice

Trademarks are the property of DENTIS CO.,LTD. or their respective owners





# Ê Cleanlant **e-Clean System**



e-Clean  
System




**DENTIS**

# DENTAL IMPLANT SOLUTION

## C o n t e n t s

- A Multi-Center Retrospective Study of DENTIS Implant System . . . . 3p
- Features of CLEANLANT . . . . . 4p
- Fixture System Selection . . . . . 5p
- Prosthetic Flow Diagrams for e-Clean . . . . . 6p
- Prosthetic Components . . . . . 8p
- Cleanlant Surgical Kit . . . . . 10p




# A Multi-Center Retrospective Study of DENTIS Implant System

Thomas K. Lee, DDS  
 Su-Kwan Kim, DDS, PhD  
 Sang-Don Joo, DDS  
 Sang-Chul Ko, DDS


**Background:** The design of DENTIS implant system, marketed since 2005 in Korea and abroad in 15 countries, is designed with the following characteristics: RBM surface treatments for time-proven osseointegration; tapered body with optimized thread designs for easy initial fixation at the time of placement surgery; 3 different abutment connection types for the same body design, allowing easier transition for the operator from existing systems in his/her armamentarium; and simplified prosthetic components. The purpose of this retrospective study was to evaluate clinical success rates for a new dental implant system called DENTIS in various private practice clinical settings.

**Methods:** 707 consecutive patients at 3 different clinical locations were treated with 1429 DENTIS implants. All 3 different abutment connection types of DENTIS Implant System (external hex connection type, octagonal conical-taper internal supra-gingival connection type, and submerged/bone level hex conical-taper internal connection type) were utilized in this study based on operator decision for each case. Implants were placed at various locations throughout the maxilla and mandible according to the treatment plan, including delayed and immediate placements after extractions. Various bone grafting procedures were done, including sinus augmentations, when clinically necessary. Patients were recalled and clinically examined at regular intervals along with radiographs to monitor clinical progression and prosthetic serviceability and stability.


## Age Distribution



## Recipients




## Anatomic Locations




**Results:** Average time since implant placements were 26 months. Average time since delivery of prosthesis was 21 months. 27 implants out of 1429 implants had to be removed before delivery of definitive restorations for various clinical failure criteria, resulting in a failure rate of 1.9%. Cumulative survival rate was 98.1%. Average age of the patient population was 52 years old at the time of implant placement surgery, while youngest patient was 16 years old and oldest patient was 87 years old. 52.7% of the patient population was female, while 47.3% was male. While maxillary molar region had the highest risk of failures anatomically, diabetes and smoking were the highest medical condition risk factors. Prosthetic complication factors such as screw loosening, cemented crowns coming-off, and porcelain fractures affected 36 implants, resulting in 4.8% prosthetic complication rate for the 26 months of this study.


## DENTIS Implant Success Rate




## Failures




## Failure Anatomic Sites



## Failure Details / Factors



## Prosthetic Complications



**Conclusion:** A new dental implant system called DENTIS is found to be performing well in various clinical practice settings in this retrospective multi-center study. Cumulative success rate of 98.1% for the DENTIS Implant System, as demonstrated in this study, compares favorably to most of the leading implant systems in the market now. Same patient population will continue be followed up in coming years for further evaluation of DENTIS implants.

# Features of CLEANLANT

e-Clean tapered

- **Fixture Mount can be used for Imperssion Coping and Free Abutment.**

Fixture mount is economical and makes it simple for dental prosthetics, Mount could be used as transfer type of impresson coping and final abutment after inserted fixture.

- **Tapered design**

Root Form Design achieved superb initial fixation strength also in places where condition of the bone is not good, makes an easy and stable surgery possible and reduces surgery time.

- **Micro-Thread for Bone Loss Prevention**


It prevents water absorption of the cortex due to bacterial infection and minimizes bone loss, breaking due to internal stress and stimulus of the cortical bone with stable initial fixation strength and wedges in thin cortex as well.

- **Safe cutting edge**

It minimizes resistance of the bone during insertion of the implant and functions to provide safe and smooth insertion.


- **RBM Surface Treatment ( $1.5 \pm 0.2 \mu\text{m}$ )**

It maximizes osseointegration.




# Fixture System Selection


e-Clean tapered



Mini		
Implant Diameter	Length	Code
Ø3.5	8.0mm	DEFM3508
	10.0mm	DEFM3510
	12.0mm	DEFM3512
	14.0mm	DEFM3514












Regular		
Implant Diameter	Length	Code
Ø4.1	8.0mm	DEFR4108
	10.0mm	DEFR4110
	12.0mm	DEFR4112
	14.0mm	DEFR4114



Wide		
Implant Diameter	Length	Code
Ø5.1	8.0mm	DEFW5108
	10.0mm	DEFW5110
	12.0mm	DEFW5112
	14.0mm	DEFW5114

# Prosthetic Flow Diagrams for e-Clean

## Cemented system

Description	Flow	Tool
Screw		 1.2 Hex Driver
<b>Abutment</b> 1) Cemented 2) Angled 3) Gold UCLA 4) Plastic		
Lab analog		
<b>Impression Coping</b> 1) Pick-up 2) Transfer		 1.2 Hex Driver
1) Cover Screw 2) Healing Abutment		 0.9 Hex Driver    1.2 Hex Driver ※Mini Only
E-Clean Fixture		

# O-Ring system





Description	Flow	Tool
1) O-Ring 2) O-Ring Retainer		
O-Ring Lab Analog		
O-Ring Abutment		<p>O-Ring Abutment Driver</p>
1) Cover Screw 2) Healing Abutment		<p>0.9 Hex Driver    1.2 Hex Driver                      ※Mini Only</p>
E-Clean Fixture		




# Prosthetic Components


## e-Clean Cemented Abutment Components

Mount	Platform Diameter	Code No.
	Mini Regular Wide	DEFMM DEFMR DEFMW


Short Mount	Platform Diameter	Code No.
	Mini Regular Wide	DEFMSM DEFMSR DEFMSW

Cover Screw	Platform Diameter	Code No.
	Mini Regular Wide	DECSM DECSR DECSW


Healing Abutment	Cuff	Code No.			
		Mini	Regular		Wide
	2mm	DEHM420	DEHR520	DEHR620	DEHW620
	3mm	DEHM430	DEHR530	DEHR630	DEHW630
	4mm	DEHM440	DEHR540	DEHR640	DEHW640
	5.5mm	DEHM455	DEHR555	DEHR655	DEHW655
	7mm	DEHM470	DEHR570	DEHR670	DEHW670



Impression Coping - Pick-Up	Type	Code No.			
		Mini	Regular		Wide
	Hex	DEIM40H	DEIR50H	DEIR60H	DEIW60H
	N-Hex	DEIM40N	DEIR50N	DEIR60N	DEIW60N

\* Impression Coping+Guide Pin Set Code : Impression Coping Code + S(ex: DEIM40HS)


Impression Coping - Transfer	Type	Code No.			
		Mini	Regular		Wide
	Hex	DETIM40H	DETIR50H	DETIR60H	DETIW60H
	N-Hex	DETIM40N	DETIR50N	DETIR60N	DETIW60N

\* Impression Coping+Guide Pin Set Code : Impression Coping Code + S(ex: DETIM40HS)


Lab Analog	Platform Diameter	Code No.
	Mini Regular Wide	DELAM DELAR DELAW

Cemented Abutment	Type	Cuff	Height(H)	Code No.			
				Mini	Regular		Wide
	Hex	6mm	1mm	DECAM416H	DECAR516H	DECAR616H	DECAW616H
			2mm	DECAM426H	DECAR526H	DECAR626H	DECAW626H
			3mm	DECAM436H	DECAR536H	DECAR636H	DECAW636H
			4mm	DECAM446H	DECAR546H	DECAR646H	DECAW646H
	N-Hex	8mm	1mm	DECAM418H	DECAR518H	DECAR618H	DECAW618H
			2mm	DECAM428H	DECAR528H	DECAR628H	DECAW628H
			3mm	DECAM438H	DECAR538H	DECAR638H	DECAW638H
			4mm	DECAM448H	DECAR548H	DECAR648H	DECAW648H
	N-Hex	6mm	1mm	DECAM416N	DECAR516N	DECAR616N	DECAW616N
			2mm	DECAM426N	DECAR526N	DECAR626N	DECAW626N
			3mm	DECAM436N	DECAR536N	DECAR636N	DECAW636N
			4mm	DECAM446N	DECAR546N	DECAR646N	DECAW646N
	N-Hex	8mm	1mm	DECAM418N	DECAR518N	DECAR618N	DECAW618N
			2mm	DECAM428N	DECAR528N	DECAR628N	DECAW628N
			3mm	DECAM438N	DECAR538N	DECAR638N	DECAW638N
			4mm	DECAM448N	DECAR548N	DECAR648N	DECAW648N


\* Abutment+Abutment Screw Set Code : Abutment Code + S(ex: DECAM416HS)

Angled Abutment	Angle	Cuff	Height(H)	Code No.			
	15°	2mm 4mm	8mm	Mini DEAAM4215 DEAAM4415	Regular DEAAR5215 DEAAR5415    DEAAR6215 DEAAR6415		Wide DEAAW6215 DEAAW6415
	25°	2mm 4mm	8mm	DEAAM4225 DEAAM4425	DEAAR5225 DEAAR5425    DEAAR6225 DEAAR6425	DEAAW6225 DEAAW6425	



\* Abutment+ Abutment Screw Set Code : Abutment Code + S(ex: DEAAM4215S)

Gold UCLA	Type	Code No.		
	Hex N-Hex	Mini DEGUMH DEGUMN	Regular DEGURH DEGURN	Wide DEGUWH DEGUWN

\* Abutment+ Abutment Screw Set Code : Abutment Code + S(ex: DEPMHS)


Plastic UCLA	Type	Code No.		
	Hex N-Hex	Mini DEPMH DEPMN	Regular DEPRH DEPRN	Wide DEPWH DEPWN


\* Abutment+ Abutment Screw Set Code : Abutment Code + S(ex: DEPMHS)


Temporary Abutment	Platform Diameter	Type	Code No.
Titanium 	Ø4.0	Titanium	DETMH
	Ø4.5	Hex	DETRH
	Ø5.5		DETWH
	Ø4.0	Titanium	DETMN
	Ø4.5	N-Hex	DETRN
	Ø5.5		DETWN
Plastic 	Ø4.0	Peek	DEPTMH
	Ø4.5	Hex	DEPTRH
	Ø5.5		DEPTWH
	Ø4.0	Peek	DEPTMN
	Ø4.5	N-Hex	DEPTRN
	Ø5.5		DEPTWN


\* Abutment+ Abutment Screw Set Code : Abutment Code + S(ex: DETMHS)

## e-Clean O-Ring Abutment Components

O-Ring Abutment	Platform Diameter	Cuff	Code No.
	Ø4.1	2mm	DEORRA20
		4mm	DEORRA40
	Ø5.1	2mm	DEORWA20
		4mm	DEORWA40

O-Ring Lab Analog	Code No.
	DOLA

O-Ring Retainer	Code No.
	DOR

O-Ring	Code No.
	ORING

# Cleanlant Surgical Kit

Probe Depth Gauge



Torque Ratchet



Code - DEK

## e - Clean tapered KIT

01 Point Drill




DPD

02 Drill Extension



DDE

Straight Twist Drill



ø 2,2

DTD2208 DTD2212  
DTD2210 DTD2214



ø 2,8


DTD2808 DTD2812  
DTD2810 DTD2814

05 Pilot Drill




ø 2,8  
DSPD28

Twist Drill




ø 3,5

DTD3508 DTD3512  
DTD3510 DTD3514



ø 4,1

DTD4108 DTD4112  
DTD4110 DTD4114



ø 5,1

DTD5108 DTD5112  
DTD5110 DTD5114

09 Parallel Pin



DPP

10 e-Clean Tapered countersink



ø 3,5 ø 4,1  
DETCS35 DETCS41  
ø 5,1  
DETCS51

11 Short Mount



Mini Regular  
DEFMSM DEFMSR  
Wide  
DEFMSW

12 Mount Driver for Machine



ø 3,5  
DMHD24  
ø 4,1  
DMHD30

13 Mount Driver for Ratchet



ø 3,5  
DRHDS24 DRHDL24  
ø 4,1  
DRHDS30 DRHDL30

14 Hex Driver



ø 0,9  
DRD09



ø 1,2

DRD12S DRD12L

## e-Clean tapered Drilling Sequence

### ø 3,5 Fixture

D4 Bone	Point drill → ø 2,2 drill → ø 2,2/2,8 pilot drill → ø 2,8 drill → ø 3,5 fixture
D2-3 Bone	Point drill → ø 2,2 drill → ø 2,2/2,8 pilot drill → ø 2,8 drill → ø 3,5 tapered drill → ø 3,5 countersink → ø 3,5 fixture
D1 Bone	Point drill → ø 2,2 drill → ø 2,2/2,8 pilot drill → ø 2,8 drill → ø 3,5 tapered drill → ø 3,5 countersink → ø 3,5 tap drill (option) → ø 3,5 fixture

### ø 4,1 Fixture

D4 Bone	Point drill → ø 2,2 drill → ø 2,2/2,8 pilot drill → ø 2,8 drill → ø 3,5 tapered drill → ø 4,1 fixture
D2-3 Bone	Point drill → ø 2,2 drill → ø 2,2/2,8 pilot drill → ø 2,8 drill → ø 3,5 tapered drill → ø 4,1 tapered drill → ø 4,1 countersink → ø 4,1 fixture
D1 Bone	Point drill → ø 2,2 drill → ø 2,2/2,8 pilot drill → ø 2,8 drill → ø 3,5 tapered drill → ø 4,1 tapered drill → ø 4,1 countersink → ø 4,1 tap drill (option) → ø 4,1 fixture

### ø 5,1 Fixture

D4 Bone	Point drill → ø 2,2 drill → ø 2,2/2,8 pilot drill → ø 2,8 drill → ø 3,5 tapered drill → ø 4,1 tapered drill → ø 5,1 fixture
D2-3 Bone	Point drill → ø 2,2 drill → ø 2,2/2,8 pilot drill → ø 2,8 drill → ø 3,5 tapered drill → ø 4,1 tapered drill → ø 5,1 tapered drill → ø 5,1 countersink → ø 5,1 fixture
D1 Bone	Point drill → ø 2,2 drill → ø 2,2/2,8 pilot drill → ø 2,8 drill → ø 3,5 tapered drill → ø 4,1 tapered drill → ø 5,1 tapered drill → ø 5,1 countersink → ø 5,1 tap drill (option) → ø 5,1 fixture

