

# DENTIS

## OneQ-SL

### Clinical Report

#### Vol.2





# DENTIS OneQ-SL Clinical Report Vol.2

## CONTENTS

1. Clinical Usefulness of DENTIS OneQ-SL Implant System in the Mandibular Posterior Region .....	04
2. Clinical Usefulness of DENTIS OneQ-SL Implant System in the Maxillary Anterior Region .....	08
3. Restoration of Posterior Tooth with DENTIS OneQ-SL Implant .....	12
4. Restoration of Posterior Tooth with DENTIS OneQ-SL Implant .....	16
5. Full Mouth Rehabilitation with Implant and Orthodontic Treatment .....	19
6. Full Mouth Rehabilitation with Implant and Orthodontic Treatment .....	23
7. Full Mouth Rehabilitation with Immediate Placement of Dental Implants .....	29
8. Restorative Driven Implant Placement with SIMPLE GUIDE .....	35
9. Implant Placement with just Two Drilling .....	40
10. Computer Guided Implant and Immediate Restoration Using SIMPLE GUIDE System .....	45
11. OneQ-SL Implantation with Lateral Sinus Window Bone Graft .....	50
12. OneQ-SL Implantation with Alveolar Ridge Splitting and GBR .....	55
13. OneQ-SL Implantation with Vertical and Horizontal Bone Augmentation .....	59
14. OneQ-SL Immediate Implantation after Tooth Extraction .....	64
15. OneQ-SL Implantation with Lateral Sinus Window Graft .....	69
16. Implant Installation and Simultaneous Horizontal Bone Graft with Hard Mesh and OneQ-SL .....	73
17. Maxillary Posterior Restoration with Sinus Elevation and Simultaneous Implant Installation .....	78
18. Sinus Elevation and Simultaneous Implant Installation with OneQ-SL .....	82

# Clinical Usefulness of DENTIS OneQ-SL Implant System in the Mandibular Posterior Region

Yun-Pyo Hong, Young-Wook Jung, Seung-Il Song, Jeong-Keun Lee  
Dept. of Oral & Maxillofacial Surgery, Ajou University School of Medicine

## Patient Information

<b>Placement Implant Area</b>			<b>Age/Sex</b>	40Y/Male		
		7				
<b>C.C.</b>	Dentigerous cyst enlarged to #37 apex #37 mobility (++)					
<b>Treatment Plan</b>	<b>1. Cyst enucleation</b> Cyst enucleation and #37, #38 extraction <b>2. Implant installation</b> - #37 implant installation after cyst enucleation in 6 months. - GBR with auto tooth bone graft material (AutoBT) <b>3. Gold crown delivery</b> #37 gold crown was delivered at 5 months after implant installation.					
<b>Material and Methods</b>	<b>Tooth No.</b>	<b>Product</b>	<b>Fixture Size</b>	<b>Torque</b>	<b>GBR</b>	<b>Bone Graft Material</b>
	#37	OneQ-SL	Ø5.2X12mm	35N	0	Auto BT

Pre-operation



Fig1. Pre-op radiography



Fig2. After cyst enucleation and #37, #38 extraction radiography




Fig3. Pre-op intraoral photo

**Intra-operation**



**Fig4.** Drilling



**Fig5.** After drilling



**Fig6.** Implant installation



**Fig7.** Healing



**Fig8.** GBR



**Fig9.** Suturing

**Table 1.** The torque and ISQ value of the implant fixture.

Placement Implant Area	Fixture Size	Torque	Immediate (buccal/lingual)	1 month (buccal/lingual)	3 months (buccal/lingual)
#37	Ø5.2X12mm	35N	68/69	58/69	69/69

The torque was measured immediately after implant 1st surgery, The ISQ values were measured immediately, 1 month and 3 months after implant 1st surgery.

### Post-operation




Fig10. Post-op radiography

### Prosthesis Delivery




Fig11. Delivery prosthesis P.A. radiograph



Fig12. Delivery prosthesis intraoral photo



Fig13. Delivery prosthesis radiography

### Conclusion

In this case, implant surgery was very successful treatment option. DENTIS OneQ implant system have good stability in unfavorable condition. The implant prosthesis show favorable state until now, after 23 months implant installation. DENTIS OneQ implant system will have quite good result so that the treatment can achieve the goal.

# Clinical Usefulness of DENTIS OneQ-SL Implant System in the Maxillary Anterior Region

Yun-Pyo Hong, Young-Wook Jung, Seung-Il Song, Jeong-Keun Lee  
Dept. of Oral & Maxillofacial Surgery, Ajou University School of Medicine

## Patient Information

<b>Placement Implant Area</b>	1   1 3	<b>Age/Sex</b>	18Y/Male			
<b>C.C.</b>	#11-23 missing due to trauma (2 years ago) Ridge atrophy on #11-23 area Ridge width : 5~6mm					
<b>Treatment Plan</b>	<b>1. GBR</b> GBR with AutoBT & Lt. ramal bone to supply insufficient width of alveolar bone. <b>2. Implant installation</b> #11, #21, #23 implant install in 1 year after bone graft					
<b>Material and Methods</b>	<b>Tooth No.</b>	<b>Product</b>	<b>Fixture Size</b>	<b>Torque</b>	<b>GBR</b>	<b>Bone Graft Material</b>
	#11	OneQ-SL	Ø4.2X10mm	17N	0	Auto BT
	#21	OneQ-SL	Ø4.2X10mm	18N	0	Auto BT
	#23	OneQ-SL	Ø4.2X10mm	17N	0	Auto BT



**Pre-operation**



Fig1. Pre-op radiography



Fig2. Pre-op intraoral photo

**Intra-operation**



**Fig3.** Drilling



**Fig4.** After drilling



**Fig5.** Implant installation



**Fig6.** Healing



**Fig7.** Suture

**Table 1.** The torque and ISQ value of the implant fixture.

Placement Implant Area	Fixture Size	Torque	Immediate (buccal/lingual)	1 month (buccal/lingual)	3 months (buccal/lingual)
#11	Ø4.2X10mm	17N	77/77	78/77	80/79
#21	Ø4.2X10mm	17N	64/63	69/69	69/69
#23	Ø4.2X10mm	17N	70/70	72/71	72/72

The torque was measured immediately after implant 1st surgery, The ISQ values were measured immediately, 1 month and 3 months after implant 1st surgery.

### Post-operation



Fig8. Post-op radiography

### Prosthesis Delivery



Fig9. Delivery prosthesis P.A. radiograph



Fig10. Delivery prosthesis intraoral photo

### Conclusion

Aesthetic restoration of anterior teeth with implant is one of the most difficult procedures to execute. Measured torque value and ISQ value for stability of DENTIS OneQ-SL and had a good result.

# Restoration of Posterior Tooth with DENTIS OneQ-SL Implant

Young-Wook Jung, Jeong-Keun Lee

Dept. of Oral & Maxillofacial Surgery, Ajou University School of Medicine

## Patient Information

<b>Placement Implant Area</b>	7	<b>Age/Sex</b>	40Y/Male		
<b>C.C.</b>	Posterior teeth missing Non smoking, non specific history of a patient visited to the office. #16 Root rest, #47 missing				
<b>Treatment Plan</b>	<b>1. #16 Bridge prosthesis delivery</b> <b>2. #47 Implant installation</b> Installed Ø5.2X10mm fixture of DENTIS OneQ-SL implant. There has no need to GBR. <b>3. Prosthesis delivery</b> After 5 months, delivered the prosthesis.				
<b>Material and Methods</b>	<b>Tooth No.</b>	<b>Product</b>	<b>Fixture Size</b>	<b>Torque</b>	<b>GBR</b>
	#47	OneQ-SL	Ø5.2X10mm	40N	X

Pre-operation



Fig1. Pre-op radiography



Fig2. Pre-op intraoral photo

**Intra-operation**



Fig3. Drilling



Fig4. Post-op intraoral photo

Table 1. The torque and ISQ value of the implant fixture.

Placement Implant Area	Fixture Size	Torque	Immediate (buccal / lingual)	1 month (buccal / lingual)	3 months (buccal / lingual)
#47	Ø5.2X10mm	40N	82/82	80/80	87/86

The torque was measured immediately after implant 1st surgery, The ISQ values were measured immediately, 1 month and 3 months after implant 1st surgery.

**Post-operation**



Fig5. Post-op radiography

## Prosthesis Delivery



Fig6. Delivery prosthesis intraoral photo



Fig7. Delivery prosthesis radiography

## Conclusion

At 3 months after the placement delivered the prosthesis, there has no specific bone absorbed at the radiography. On 2-year follow-up, there has no abnormal finding and maintain with stable condition. DENTIS OneQ-SL Implant has a taper-straight design for high initial stability and surface is S.L.A surface for faster osseointegration. S.L.A surface is already known as an excellent surface through experiments.

In this study also indicate high ISQ and Insertion torque value. And until now, there has no abnormal finding and maintain with stable condition.

But within the limitation of this study, only installed at the posterior mandible site. Therefore, I offer more variety clinical study in the future.

# Restoration of Posterior Tooth with DENTIS OneQ-SL Implant

Young-Wook Jung, Jeong-Keun Lee

Dept. of Oral & Maxillofacial Surgery, Ajou University School of Medicine

## Patient Information

<b>Placement Implant Area</b>	7	<b>Age/Sex</b>	60Y/Female			
<b>C.C.</b>	Posterior teeth missing Non smoking, non specific history of a patient visited to the office. #47 missing					
<b>Treatment Plan</b>	<p><b>1. #47 GBR</b> GBR with AutoBT &amp; Bovine bone to supply insufficient width of alveolar bone. And covered with a resorbable collagen membrane</p> <p><b>2. #47 Implant installation</b> Installed Ø5.2X10mm fixture of DENTIS OneQ-SL implant after 5 months. There has no additional GBR.</p> <p><b>3. Prosthesis delivery</b> After 5 months, delivered the prosthesis.</p>					
<b>Material and Methods</b>	<b>Tooth No.</b>	<b>Product</b>	<b>Fixture Size</b>	<b>Torque</b>	<b>GBR</b>	<b>Bone Graft Material</b>
	#47	OneQ-SL	Ø5.2X10mm	35N	0	AutoBT, Bovine bone, Resorbable collagen membrane



**Pre-operation**




Fig1. Pre-op radiography

**Intra-operation**




Fig2. GBR

Table 1. The torque and ISQ value of the implant fixture.

Placement Implant Area	Fixture Size	Torque	Immediate (buccal/lingual)	1 month (buccal/lingual)	3 months (buccal/lingual)
#47	Ø5.2X10mm	35N	67/67	72/75	75/79

The torque was measured immediately after implant 1st surgery, The ISQ values were measured immediately, 1 month and 3 months after implant 1st surgery.

**Post-operation**




Fig3. Post-op radiography

### Prosthesis Delivery



Fig4. Post-op intraoral photo

### Conclusion

At 5 months after the implant placement delivered the prosthesis, there has no specific bone absorbed at the radiography after 2 months. On 2-year follow-up, there has no abnormal finding and maintain with stable condition. In this study also indicate high ISQ and Insertion torque value. And until now, there has no abnormal finding and maintain with stable condition.

## Full Mouth Rehabilitation with Implant and Orthodontic Treatment

Young-Lin Cho<sup>1</sup>, U-Ju Yang<sup>2</sup>, Hye-Ju Kang<sup>3</sup>, Jong-Chul Park<sup>4</sup>

<sup>1</sup> Ph.D., Director, Withus Dental Clinic

<sup>2</sup> Director, Yonsei Bon Dental Clinic

<sup>3</sup> M.S.D. Director, Shinwoo Dental Clinic

<sup>4</sup> Director, Withus Dental Clinic

### Patient Information

<b>Placement Implant Area</b>	<table border="1"> <tr> <td>6 4</td> <td>5 6</td> </tr> <tr> <td>7</td> <td>7</td> </tr> </table>	6 4	5 6	7	7	<b>Age/Sex</b>	49Y/Female
6 4	5 6						
7	7						
<b>C.C.</b>	Tooth space due to missing teeth, Missing teeth of #16, #26, #37, Root rest of #16, #47						
<b>Treatment Plan</b>	<p>Patient had spaces on the anterior area and wanted overall treatment. And she also complained loss of masticatory function because of several posterior missing teeth. We made the treatment plan for masticatory function recovery, occlusion, and appearance improvement.</p> <p>1 year after Implant Installation, orthodontic treatment was finished.</p> <p>Delivery #47 zirconia crown, #16 x 14 and #25, 26, 37 PFM crown and bridge</p>						
<b>Material and Methods</b>	<b>Tooth No.</b>	<b>Product</b>	<b>Fixture Size</b>	<b>Torque</b>	<b>Significant Content</b>		
	#37	OneQ-SL	Ø3.7X10mm	50N	-		
	#14	OneQ-SL	Ø3.7X8mm	50N	Immediate Placement		
	#47	OneQ-SL	Ø4.2X8mm	50N	Immediate Placement, BAOSFE		
	#16	OneQ-SL	Ø5.7X8mm	50N	Immediate Placement, BAOSFE		
	#25	OneQ-SL	Ø7.0X8mm	50N	Immediate Placement		
	#26	OneQ-SL	Ø4.2X8mm	50N	Lateral sinus elevation, Submerging		

**Pre-operation**



Fig1. Pre-op extraoral photo



Fig2. Pre-op radiography



Fig3. Pre-op intraoral photo

---

### Bracketing and Implant Installation



Fig4. 1.5 months after orthodontic treatment

### Post-operation




Fig5. Post-op radiography

### Prosthesis Delivery



Fig6. Final prosthesis delivery extraoral photo



Fig7. Final prosthesis delivery intraoral photo



Fig8. Final prosthesis delivery radiography

### Conclusion

Patient was treated with immediate implant accompanied by sinus augmentation and orthodontic treatment. An implant of #16 was placed with crestal sinus lift, an Implant of #25 was immediately placed and accompanied by sinus lift, and an implant of #26 was placed with lateral sinus lift technic. If all procedures were successful, there are no prognosis differences among those surgical techniques. At the first panoramic x-ray, we can see the thickening of sinus membrane. But after the removing the tooth causing trouble and placing the implant, the sinus healed without other treatment.

## Full Mouth Rehabilitation with Implant and Orthodontic Treatment

Young-Lin Cho<sup>1</sup>, U-Ju Yang<sup>2</sup>, Hye-Ju Kang<sup>3</sup>, Jong-Chul Park<sup>4</sup>

<sup>1</sup> Ph.D., Director, Withus Dental Clinic

<sup>2</sup> Director, Yonsei Bon Dental Clinic

<sup>3</sup> M.S.D. Director, Shinwoo Dental Clinic

<sup>4</sup> Director, Withus Dental Clinic

### Patient Information

<b>Placement Implant Area</b>	7 6 7 6 2	5 6 7 4 5 7	<b>Age/Sex</b>	49Y/Female	
<b>C.C.</b>	#33-43 Bridge, Missing teeth of #32-43, 36, 45, Crowding				
<b>Treatment Plan</b>	Patient complained about bimaxillary protrusion, loss of masticatory function and crowding. 6 months after Implant Installation, delivery #17-16 PFM Bridge, #25-27 PFM Bridge. 11 months after Implant Installation, orthodontic treatment was finished. Delivery #4746 x zirconia bridge, #43 42 x x x 34 PFM bridge, #35 x 37 zirconia bridge				
<b>Material and Methods</b>	<b>Tooth No.</b>	<b>Product</b>	<b>Fixture Size</b>	<b>Torque</b>	<b>Significant Content</b>
	#42	OneQ-SL	Ø2.0X10mm	50N	Immediate Placement
	#34	OneQ-SL	Ø4.2X10mm	50N	Immediate Placement
	#35	OneQ-SL	Ø4.7X8mm	50N	-
	#37	OneQ-SL	Ø3.7X8mm	50N	-
	#46	OneQ-SL	Ø4.2X8mm	50N	-
	#47	OneQ-SL	Ø4.2X8mm	50N	-
	#17	OneQ-SL	Ø4.7X8mm	50N	OSFE, Submerging
	#16	OneQ-SL	Ø4.2X8mm	50N	OSFE, Submerging
	#25	OneQ-SL	Ø3.7X8mm	50N	Immediate Placement
	#26	OneQ-SL	Ø4.2X8mm	50N	BAOSFE, Submerging
#27	OneQ-SL	Ø3.7X8mm	50N	BAOSFE, Submerging	

**Pre-operation**



Fig1. Pre-op extraoral photo



Fig2. Pre-op radiography

**Mandibular Post-operation**




Fig3. Mandibular post-op radiography



## Mandibular Prosthesis Delivery



Fig4. Mandibular temporary prosthesis delivery radiography



Fig5. Mandibular anterior temporary prosthesis delivery intraoral photo



Fig6. Mandibular final prosthesis delivery radiography



Fig7. Mandibular final prosthesis delivery intraoral photo

### Maxillary Post-operation



Fig8. Maxillary post-op radiography

### Maxillary Prosthesis Fabrication



Fig9. Maxillary final prosthesis fabrication

**Prosthesis Delivery**



**Fig10.** Final prosthesis delivery intraoral photo



**Fig11.** Final prosthesis delivery radiography



**Fig12.** Final prosthesis delivery extraoral photo

## Conclusion

Patient had overall periodontal disease and especially the lower jaw has decreased alveolar ridge and limited gingiva by edentulous state for a long time. Upper jaw has a very poor periodontal condition and needs full mouth rehabilitation. After upper posterior teeth extraction, During the orthodontic treatment, the lower teeth have restored. upper posterior area's first treatment plan was to implant accompanied by lateral sinus lift , but the plan was altered to crestal sinus lift with GBR because of improper operation view and difficulty to approach. In like this case, if proper treatment is performed, we look forward to the good prognosis through crestal approach procedure. And her appearance has improved and had stable occlusion.

## Full Mouth Rehabilitation with Immediate Placement of Dental Implants

Young-Lin Cho<sup>1</sup>, U-Ju Yang<sup>2</sup>, Hye-Ju Kang<sup>3</sup>, Jong-Chul Park<sup>4</sup>

<sup>1</sup> Ph.D., Director, Withus Dental Clinic

<sup>2</sup> Director, Yonsei Bon Dental Clinic

<sup>3</sup> M.S.D. Director, Shinwoo Dental Clinic

<sup>4</sup> Director, Withus Dental Clinic

### Patient Information

<b>Placement Implant Area</b>	6 4 3 7 5 2	3 5 6 3 5 6	<b>Age/Sex</b>	72Y/Male	
<b>C.C.</b>	Root rest of #12-17, 22, 23, 27, 34-37, 45, Missing teeth of #24-26, 47				
<b>Treatment Plan</b>	<p>A patient visited to consult about overall dental treatment. The patient's bite collapsed and periodontal condition was very poor. This case is treated with immediate placement of implant.</p> <p>#42 x x x 33 x 35 36 temporary bridge, 1 month after implant installation          #42 x x x 33 x 35 36 PFM bridge delivery, 3 months after implant installation          #16 x 14 zirconia bridge, 3 months after implant installation          #47 x 45 zirconia bridge, 3 months after implant installation          #13 x x x x 23 x 25 26 temporary bridge, 4 months after implant installation          #13 x x x x 23 x 25 26 zirconia bridge, 6 months after implant installation</p>				
<b>Material and Methods</b>	<b>Tooth No.</b>	<b>Product</b>	<b>Fixture Size</b>	<b>Torque</b>	<b>Significant Content</b>
	#42	OneQ-SL	Ø4.7X10mm	50N	Immediate Placement
	#33	OneQ-SL	Ø4.2X12mm	50N	Immediate Placement
	#35	OneQ-SL	Ø4.7X10mm	50N	Immediate Placement
	#36	OneQ-SL	Ø4.2X10mm	50N	Immediate Placement
	#25	OneQ-SL	Ø4.2X10mm	50N	Immediate Placement
	#26	OneQ-SL	Ø4.7X7mm	50N	OSFE, submerging
	#47	OneQ-SL	Ø4.7X10mm	50N	Immediate Placement
	#45	OneQ-SL	Ø4.2X10mm	50N	Immediate Placement
	#16	OneQ-SL	Ø4.2X8mm	50N	Immediate Placement
	#14	OneQ-SL	Ø4.2X12mm	50N	Immediate Placement
	#13	OneQ-SL	Ø4.2X12mm	50N	Immediate Placement
	#23	OneQ-SL	Ø4.2X12mm	50N	Immediate Placement, Submerging

**Pre-operation**



Fig1. Pre-op extraoral photo




Fig2. Pre-op radiography



Fig3. Pre-op intraoral photo

Intra-operation



Fig4. Mandibular incisor implant installation



Fig5. Mandibular left side immediately implant installation after extraction




Fig6. #45, #47 root rest extraction



Fig7. #45, #47 implant installation



Fig8. #45 GBR



Fig9. #45 membrane



**Fig10.** #45, #47 Suture



**Fig11.** #26 Extraction



**Fig12.** #25, #26 Implant Installation



**Fig13.** #25, #26 Suture

**Post-operation**



**Fig14.** Upper right side immediately implant installation after extraction



**Fig15.** Post-op radiography



**Fig16.** Healing period





Prosthesis Delivery



Fig17. Provisional temporary prosthesis delivery



Fig18. Final prosthesis delivery



Fig19. Final prosthesis delivery extraoral photo



Fig20. Final prosthesis delivery radiography

## Conclusion

A patient visited the office to consult about overall dental treatment. This patient's overall bite collapsed and periodontal condition was also very poor. It'll be very long time consuming to implant into healed ridge after extraction teeth. So, with the judgement that immediate implant after extraction causes no trouble and reduces the treatment period remarkably, we made a plan for full mouth rehabilitation using immediate placement and loading of implant. After treatment, the patient had generally harmonious jaw relationship.

As a result, we reduced the treatment period remarkably and minimized the loss of dental esthetics and function for the patient. Consequently, dentist and patient were satisfied with treatment result. In this way, immediate placement and loading of implant accompanied by properly performed bone graft is thought to be worth treatment option which is fully predictable.

Still, There are many different opinions about immediate placement of implant into fresh extraction site. But based from my 15-year-long experience of immediate implant surgery after extraction, It has showed good result without big problem. Of course, I'm not suggesting applying it to all cases but selectively applied, it is fully predictable and worth of choosing treatment option.

---

# Restorative Driven Implant Placement with SIMPLE GUIDE

Marco Tulio Alzaga Vega  
Director GDIA Mexico

## Patient Information

<b>Placement Implant Area</b>	<table border="1"> <tr> <td>x x x</td> <td>x x</td> </tr> </table>		x x x	x x	<b>Age/Sex</b>	50Y/Female			
x x x	x x								
<b>C.C.</b>	The Patient found it uncomfortable to eat by her lower denture. Excessive pressure was placed in the area of the mental foramen. Radiographs reflected loss of the majority of her bone structure.								
<b>Treatment Plan</b>	<p>It was decided to place 5 implants in the area between both metal foramens to secure for a mandibular implant supported screw retained prosthesis. During this guided surgical process, we followed the "SIMPLE GUIDE Protocol".</p> <p><b>1. Scanning Stage:</b></p> <ol style="list-style-type: none"> <li>1) Fabricate Scan Appliance (Using wax)</li> <li>2) CT scan the patient while wearing the scan appliance.</li> <li>3) Model scanning for the model with the scan appliance.</li> <li>4) Model scanning for the model without the scan appliance.</li> </ol> <p><b>2. Planning Stage:</b></p> <ol style="list-style-type: none"> <li>1) Loading 3 of CT and model scanning data on the BlueSkyPlan</li> <li>2) Combing 3 of data</li> <li>3) Designing</li> </ol> <p><b>3. Printing Stage:</b></p> <ol style="list-style-type: none"> <li>1) Printing Guide Stent</li> </ol> <p><b>4. Surgical Operation Stage:</b></p> <ol style="list-style-type: none"> <li>1) Wearing Guide Stent</li> <li>2) Drilling with SIMPLE GUIDE Kits</li> <li>3) Removed Guide Stent</li> <li>4) Drilling with OneQ Kits</li> <li>5) Implant Installation</li> </ol>								
<b>Material and Methods</b>	<b>Tooth No.</b>	<b>Product</b>	<b>Fixture Size</b>	<b>Torque</b>	<b>Guide System</b>	<b>GBR</b>	<b>Bone Graft Material</b>		
	Mandibular All-on-five	OneQ-SL	Ø3.7X10mm	45N	SIMPLE GUIDE	0	Allograft, Berkeley, Collagen membrane, Rapiderm		

**Pre-operation**



**Fig1.** Pre-op extraoral photo



**Fig2.** Pre-op radiography



**Fig3.** Pre-op intraoral photo



**Fig4.** Try-in of provisional restoration

### Guide Stent Fabrication Procedure



Fig5. Impressions are taken and a master model is obtained



Fig6. A diagnostic wax-up created with 6-8 fiducial markers in place.






Fig7. CBCT of the patient with the scan appliance



Fig8. Optical scan of model without scan appliance.




Fig9. Optical scan of model with scan appliance.




Fig10. Digital planning (BlueSky Software) combines CT imaging and model scans.




Fig11. Digital image of the guide



Fig12. 3D printed guide stent with open metallic sleeves.



Fig13. Guide stent in place

### Intra-operation



**Fig14.** Mid-crestal incision



**Fig15.** Using a piezo to clean the ridge after extractions.



**Fig16.** Drilling of the anchor screw.



**Fig17.** Anchor screws stabilize the guide into place.



**Fig18.** Drilling sequence with the guide in place.



**Fig19.** Implants in place obtaining over 45Ncm Torque



**Fig20.** OneQ-SL Ø3.7 X 10mm were place on all five osteotomy

### Conclusion

Choose the system with the appropriate surgical and prosthetic options necessary for each case. Various systems of guided surgery are available, all producing favorable results. The SIMPLE GUIDE protocol is an efficient and predictable method, providing the operator with the option to select the most appropriate implant system necessary for specific surgical or prosthetic concerns.

# Implant Placement with just Two Drilling

Hyeon-Min Kim <sup>1</sup>, Jae-Young Lyu <sup>1</sup>, Jin-Young Cho <sup>1</sup>, Sung-Won Yang <sup>2</sup>, Woo-Yul Lee <sup>2</sup>

<sup>1</sup> Prof., Dept. OMFS, Gil Medical Center, Gachon University

<sup>2</sup> DDS, Dept. OMFS, Gil Medical Center, Gachon University

## Patient Information

<b>Placement Implant Area</b>	6		<b>Age/Sex</b>	45Y/Male		
<b>C.C.</b>	The patient feels discomfort at chewing. And he complains about caries of #17, mobility of #46 tooth.					
<b>Treatment Plan</b>	<p><b>1. Preparation stage for implant installation</b> 1) Ext. of #46</p> <p><b>2. Surgical planning stage</b> 1) Analysis of bone width &amp; height in CBCT 2) Decision of implant diameter, length, path</p> <p><b>3. Operation stage</b> 1) 1<sup>st</sup> surgery : implant installation(#46) 2) 2<sup>nd</sup> surgery : healing abutment(#46)</p> <p><b>4. Prosthetic stage</b> 1) Impression : transfer imp. (#46) 2) Coping try-in(#46) 3) Temp. setting(#46) 4) Final setting(#46)</p>					
<b>Material and Methods</b>	<b>Tooth No.</b>	<b>Product</b>	<b>Fixture Size</b>	<b>Torque</b>	<b>1<sup>st</sup> OP ISQ</b>	<b>2<sup>nd</sup> OP ISQ</b>
	#46	OneQ-SL	Ø4.7X10mm	>30 Ncm	75	76



Pre-operation



**Fig1.** Pre-op intraoral photo. Gingival healing state after #46 ext. was identified.



**Fig2.** Pre-op radiography. Bony healing state after #46 ext. is identified.



**Fig3.** Diameter, length and path of implant were considered using 3D implant simulation.

### Implant 1<sup>st</sup> surgery



**Fig4.** Crestal and distal vertical incision were performed. And the flap was elevated.



**Fig5.** Just two drilling (initial and final drill) was done. #46 implant installation was performed.



**Fig6.** Primary stability could be abstained.



**Fig7.** Implant installation

### Post- Implant 1<sup>st</sup> surgery



**Fig8.** Diameter, length and path of implant were considered using 3D implant simulation.

**Implant 2<sup>nd</sup> surgery**



**Fig9.** Pre-2<sup>nd</sup> op intraoral photo. Good gingival healing after 1st op.



**Fig10.** Implant 2<sup>nd</sup> surgery on #46 area after crestal incision(#46 ext. site) and flap elevation were done.




**Fig11.** Healing abutment was tightened.

**Post- Implant 2<sup>nd</sup> surgery**



**Fig12.** Post-2<sup>nd</sup> op (#46) radiography



**Fig13.** Healing abutment was tightened on #46 fixture.

### Prosthesis Delivery




Fig14. Impression taking was done after 2weeks (#46)



Fig15. Coping(#46) Check



Fig16. Final Restoration (#46)

### Conclusion

The aim of this case report was the evaluation of usefulness of implant stability after skipped drill order. All the procedure was same with common implant procedure except 1<sup>st</sup> surgery. The results showed good implant stability at both 1<sup>st</sup> and 2<sup>nd</sup> surgery stage.

# Computer Guided Implant and Immediate Restoration Using SIMPLE GUIDE System

Jin-yong Cho

Clinical assistant professor, Department of oral and maxillofacial surgery, Gachon University Gil Medical Center

## Patient Information

<b>Placement Implant Area</b>	3	<b>Age/Sex</b>	50Y/Female
<b>C.C.</b>	Prolonged retained primary tooth and its mobility		
<b>Treatment Plan</b>	In this report, we would present a patient that was treated by implant and immediate restoration using SIMPLE GUIDE System, i.e. simple and minimal-cost computer guided implant system.		
<b>Material and Methods</b>	<ol style="list-style-type: none"> <li>1. #53 primary tooth was dropped out spontaneously</li> <li>2. Preoperative scan data and simulation of #53 extraction using Meshmixer software</li> <li>3. Planning and guide design using Blue sky plan III software.</li> <li>4. Surgical guide printed with ZENITH 3D printer and temporary crown fabrication in the dental model which was made using surgical guide.</li> <li>5. Application of surgical guide and implant installation with OneQ-SL Ø3.7 X 10mm</li> <li>6. Guided system, SIMPLE GUIDE Kits</li> </ol>		

**Pre-operation**



**Fig1.** Pre-op radiography. It showed #13 ectopic impaction and prolonged retained primary tooth on #53. Narrowing inter-root distance between #12 and #14 was noted.



**Fig2.** Pre-op intraoral photo. #53 primary tooth was dropped out spontaneously between initial visit and operation day.

### Guide Stent Fabrication Procedure



Fig3. Preoperative scan data (Lt) and simulation of #53 extraction using Meshmixer software. (Rt)





Fig4. Planning and guide design using Blue sky plan III software. A. Decide the implant position corresponding ideal prosthesis. B Merge the CT and intraoral scan data. C. Design the guide outline. D. Automatically fabricated surgical guide.



Fig5. Surgical guide printed with ZENITH 3D printer and temporary crown fabrication in the dental model which was made using surgical guide.

**Intra-operation**




**Fig6.** Application of surgical guide and implant installation.

**Prosthesis Delivery**



**Fig7.** Immediate temporary crown is restored with favorable suitability.



**Fig8.** Difference between planning (red line) and actual installed implant. Coronal view (Lt) and axial view (Rt).



## Case 10

---

**Table 1.** Accuracy of guided implant. Coordinate the planning implant and actual installed implant at the center of top and bottom. Calculate the difference between two coordinates.

	<b>Top</b>	<b>Bottom</b>
d-x	0.49	0.42
d-y	0.81	0.38
d-z	0.65	0.13

mm

	<b>Axis difference</b>
3D	1.79
2D	1.8

degree

### Conclusion

In this case, implant and immediate restoration in the maxillary anterior zone using SIMPLE GUIDE system is considered as simple and reliable method with high accuracy and low-cost.

# OneQ-SL Implantation with Lateral Sinus Window Bone Graft

Jong-Hun Jeong  
Light Dental Office, Gyeongju, Republic of Korea

## Patient Information

<b>Placement Implant Area</b>	7		<b>Sex</b>	Female		
<b>C.C.</b>	Pain on upper left molar area					
<b>P.I.</b>	1. #26 and #27 old bridge 2. Peri-implant radiolucency seen at tooth #27.					
<b>Treatment Plan</b>	1. Removal of #26 and #27 old bridge 2. Extraction of #27 3. Extract #27 and waiting 2 months for delayed implantation. 4. Implantation using GBR via a lateral sinus window with vertical bone augmentation.					
<b>Material and Methods</b>	1. The lateral sinus window is opened. 2. Bone grafting within the sinus with Xeno-bone. 3. Placement completed, using an OneQ-SL Ø4.7 X 10mm implant. 4. Perform vertical bone augmentation with autogenous bone in the gap between the implant fixture and the remaining extraction socket. 5. Schedule secondary patient visit. 6. A final restoration is placed after loading the provisional restoration.					
	<b>Tooth No.</b>	<b>Product</b>	<b>Fixture Size</b>	<b>Torque</b>	<b>Sinus OP</b>	<b>Bone Graft Material</b>
	#27	OneQ-SL	Ø4.7X10mm	30 Ncm	0	Xeno-bone

Pre-operation



Fig1. Pre-op #27 P.A. radiography



Fig2. Post extraction P.A. radiography




Fig2. Pre-op radiography



Fig3. Pre-op intraoral photo

---

**Implant 1<sup>st</sup> surgery**



**Fig5.** Lateral sinus window opening



**Fig6.** Sinus floor elevation



**Fig7.** Sinus bone grafted



**Fig8.** Sinus window closed with lateral bone wall



**Fig9.** #27 implantation : OneQ-SL Ø4.7x10mm (11Ncm)



**Fig10.** Bone augmentation with autogenous bone

**Post- Implant 1<sup>st</sup> surgery**



**Fig11.** Post-op radiography



**Fig12.** Post-op P.A. radiography after 2 weeks



**Fig13.** Post-op intraoral photo after 2 weeks



**Fig14.** Post-op P.A. radiography after 4 months



**Fig15.** Post-op intraoral photo after 4 months

### Implant 2<sup>nd</sup> surgery



Fig16. 2<sup>nd</sup> op intraoral photo



Fig17. Healing abutment setting after 2<sup>nd</sup> op

### Post- Implant 2<sup>nd</sup> surgery




Fig18. P.A. radiography after setting a healing abutment



Fig19. Abutment setting and loading P.A. radiography after 22 weeks

### Conclusion

In this situation, we sought to achieve early stability using under drilling techniques. It was expected that the OneQ-SL implant fixture's thread design and great surface treatment of SLA would allow for bone augmentation into the sinus and implantation of the implants to be performed simultaneously. OneQ-SL implant fixture's strong thread, tapered design allowed it to fit securely into the bone, achieving an early fixed torque. Optimal osseointegration was seen during the second appointment, 18 weeks later. 22 weeks after the surgery, 30Ncm of force was used when connecting the abutment and loading of the provisional restoration. No challenges were found when applying enough torque force. Radiographically, osseointegration condition between the thread of the implant fixture looked optimal, along with satisfactory marginal bone conditions. Overall, the SLA surface's osseointegration capabilities were satisfactory.

# OneQ-SL Implantation with Alveolar Ridge Splitting and GBR

Jong-Hun Jeong  
 Light Dental Office, Gyeongju, Republic of Korea

## Patient Information

<b>Placement Implant Area</b>			<b>Sex</b>	Female
		6 7		
<b>C.C.</b>	Prothodontic rehabilitation on lower left molar area			
<b>P.I.</b>	1. Vertical eruption of #26 and #27 due to missed #36 and #37 2. Narrow alveolar ridge on #36 and #37 area			
<b>Treatment Plan</b>	1. #26, #27: 2 unit bridge after Endo Tx. 2. Ridge splitting and GBR – resorbable membrane 3. Submerged implantation 4. Need to 2nd OP after 2-4 months on implantation (Consider APF flap if need) 5. Final restoration after loading of provisional restoration			
<b>Material and Methods</b>	1. Submerged implantation on #36 and #37 area with OneQ-SL Ø4.2x10mm 2. Ridge splitting and GBR with resorbable membrane			

**Pre-operation**



**Fig1.** Pre-op radiography



**Fig2.** Pre-op intraoral photo

**Implant 1<sup>st</sup> surgery**



**Fig3.** Alveolar ridge splitting (Bone sawing and initial hole formation)



**Fig4.** Alveolar ridge bone splitting and spreading



**Fig5.** Hole formation for implantation




**Fig6.** #36 Implantation OneQ-SL Ø4.2x10mm (10Ncm)



**Post- Implant 1<sup>st</sup> surgery**



**Fig7.** Post-op radiography



**Fig8.** Post-op P.A. radiography



**Fig9.** Post-op 11 weeks ridge foam

**Implant 2<sup>nd</sup> surgery**



**Fig10.** 2<sup>nd</sup> op intraoral photo



**Fig11.** Suture

### Post- Implant 2<sup>nd</sup> surgery



Fig12. 2<sup>nd</sup> op P.A. radiography after op 11 weeks

### Prosthesis Delivery



Fig13. Abutment setting (Provisional restoration / Loading) P.A. radiography after op 14 weeks



Fig14. Final restoration P.A. radiography after op 19 weeks


### Conclusion

Here, implantation was performed using alveolar ridge splitting and GBR. During installation of implants, initial stability is key, especially when ridge splitting occurs during immediate implantation following an extraction. On the radiograph, osseointegration condition between the thread of the implant fixture were optimal, with satisfactory marginal bone conditions. Overall, SLA surfaces' osseointegration capabilities are extremely satisfactory, even with initial ridge splitting occurring.

## OneQ-SL Implantation with Vertical and Horizontal Bone Augmentation

Jong-Hun Jeong  
Light Dental Office, Gyeongju, Republic of Korea

### Patient Information

<b>Placement Implant Area</b>		<b>Sex</b>	Female
<b>C.C.</b>	Pain on upper left molar area		
<b>P.I.</b>	A periodontal lesion is noted radiographically surrounding the root of tooth #27, which is found to be in a floating state.		
<b>Treatment Plan</b>	<ol style="list-style-type: none"> <li>1. Recommended removal of bridge and extraction of -#25 and #27.</li> <li>2. Delayed implantation and vertical bone augmentation (GBR) in area of #27.</li> <li>3. Submerged implantation at #26 and #27.</li> <li>4. Schedule patient for secondary visit.</li> </ol>		
<b>Material and Methods</b>	<ol style="list-style-type: none"> <li>1. Wait 2-3 months following extraction of #27.</li> <li>2. Delayed submerged implantation at #26 and #27.</li> <li>3. Vertical augmentation (GBR) with autogenous bone and non-resorbable membrane</li> <li>4. 2<sup>nd</sup> OP scheduled in 4 months.</li> </ol>		

Pre-operation



**Fig1.** Pre-op #26-27 P.A. radiography. Observation widely bone loss #27 periimplant.



**Fig2.** Postextraction of #27 10 weeks. Observed that #27 extraction socket remain widely. And need a vertical bone augmentation


**Implant 1<sup>st</sup> surgery**



**Fig3.** After extraction 10 weeks. Observation soft tissue healed well



**Fig4.** . Wide bone defect on #27 extraction socket area




**Fig5.** Hole formation for implantation (under drilling and bone condensation with osteotome)



**Fig6.** Implant placement #26 OneQ-SL Ø4.7x12 (55Ncm), #27 OneQ-SL Ø4.7x12 (42Ncm)




**Fig7.** Bone defect at buccal observation



**Fig8.** Bone defect at palatal observation



**Fig9.** GBR




**Fig10.** Non-resorbable membrane coverage

**Post- Implant 1<sup>st</sup> surgery**



**Fig11.** Post-op radiography



**Fig12.** Post-op P.A. radiography



**Fig13.** Post-op after 10 weeks

**Implant 2<sup>nd</sup> surgery**




**Fig15.** Bone filled between the implant fixture thread




**Fig16.** Observed bone formation well

### Post- Implant 2<sup>nd</sup> surgery



**Fig14.** Post-op P.A. radiography after 2<sup>nd</sup> op 14 weeks

### Prosthesis Delivery



**Fig17.** Post-op 5 months (Loading). Setting Abutment (30Ncm). Loading by provisional restoration



**Fig18.** Final restoration P.A. radiography



**Fig19.** Final restoration intraoral photo


### Conclusion

With enough time for osseointegration given for the area of bone augmentation, there was concern that bone quality might be reduced. A provisional restoration was loaded with adequate torque achieved during connection. Radiographically, osseointegration conditions between the implant thread and bone (along with marginal bone conditions) appeared to be satisfactory.

# OneQ-SL Immediate Implantation after Tooth Extraction

Jong-Hun Jeong  
 Light Dental Office, Gyeongju, Republic of Korea

## Patient Information

<b>Placement Implant Area</b>		<b>Sex</b>	Female
<b>C.C.</b>	Pain on lower left molar area		
<b>P.I.</b>	1. #34-36 distal cantilever bridge (#36 pontic) 2. #34 – horizontal tooth fracture #35 – severe dental caries and root rest #36 – missing		
<b>Treatment Plan</b>	1. Immediate implantation after #34 and #35 extraction and implantation on #36 area 2. Bone augmentation (GBR) on gap between implant fixture and extraction socket 3. Non-submerged OP with healing abutment connection		
<b>Material and Methods</b>	1. Extraction of #34 and #35 2. Immediate implantation on #34-36 area 3. Bone augmentation (GBR) was performed using autogenous bone in the area between the implant fixture and residual extraction socket. 4. Healing abutment connection		



Pre-operation



Fig1. Pre-op radiography



Fig2. #34-36 Distal cantilever bridge



Fig3. #34 Horizontal Fx.





Fig2. #35 Severe dental caries (Root rest)

**Intra-operation**



**Fig5.** Extraction socket on #34-35



**Fig6.** Hole formation for immediate implantation



**Fig7.** #34 Ø4.2x12mm [21Ncm], #35 Ø4.2x12mm [25Ncm], #36 Ø4.2x10mm [27Ncm]



**Fig8.** Bone augmentation (GBR). Autogenous bone placed in space between implant fixture and extraction socket.



**Fig9.** GBR (unused membrane)



**Fig10.** Suture

Post-operation



Fig11. Post-op radiography




Fig12. Post-op P.A. radiography

## Prosthesis Delivery




Fig13. Post-op 3 months after abutment connection (30Ncm)




Fig14. Final restoration P.A. radiography



Fig15. Final prosthesis delivery after 14 weeks intraoral photo

## Conclusion

In this case, a provisional restoration was loaded first, with adequate torque force on the connecting abutment. Radiographically, osseointegration conditions between the implant's SLA-treated threading surfaces and marginal bone were satisfactory .

## OneQ-SL Implantation with Lateral Sinus Window Graft

Jong-Hun Jeong  
 Light Dental Office, Gyeongju, Republic of Korea

### Patient Information

<b>Placement Implant Area</b>	5 6	<b>Sex</b>	Female
<b>C.C.</b>	Prothodontic rehabilitation on upper molar missed dentition		
<b>P.I.</b>	1. Missed tooth of #16, #17, #25, #26 2. Floaying tooth of #14 and #15 3. Severe rhematization of both sinus floors.		
<b>Treatment Plan</b>	1. Extraction of #14 and #15 – delayed implantation with GBR 2. Implantation at #15 and #16 was accessed via lateral sinus window grafting with Xeno-bone. 3. Delayed implantation at #14 and #16 via lateral sinus window grafting with Xeno-bone and vertical bone augmentation using Allobone and membrane.		
<b>Material and Methods</b>	<b>Prosedure on #25 and #26 area</b> 1. Lateral sinus window opening 2. Sinus bone grafting with Xeno-bone 3. Final restoration after provisional restoration		

Pre-operation




Fig1. Pre-op radiography

**Intra-operation**



**Fig2.** Lateral sinus window opening



**Fig3.** Closed sinus window with bone wall after sinus bone graft



**Fig4.** Hole formation for implantation



**Fig5.** OneQ-SL implant implantation, #25 Ø4.2x10mm (21Ncm), #26 Ø4.2x10mm (21Ncm)



**Fig6.** Suture

**Post-operation**



Fig7. Post-op radiography




Fig8. Post-op P.A. radiography after 5 months

**Prosthesis Delivery**



Fig9. Connected abutment P.A. radiography



Fig10. Temporary Crown setting





# Implant Installation and Simultaneous Horizontal Bone Graft with Hard Mesh and OneQ-SL

Dae-Sung Kim  
 DDS, MSD.Ph.D, Director, MediPlant D.C/MEDI Academy of Implant Dentistry

## Patient Information

<b>Placement Implant Area</b>	7 6	<b>Age/Sex</b>	48Y/Female		
<b>C.C.</b>	Long term edentulous posterior region				
<b>Treatment Plan</b>	Implant installation and simultaneous horizontal bone graft with hard Mesh				
<b>Material and Methods</b>	<b>Tooth No.</b>	<b>Product</b>	<b>Fixture Size</b>	<b>GBR</b>	<b>Bone Graft Material</b>
	#46	OneQ-SL	Ø4.2X12mm	0	Cerasorb M Titanium Mesh
	#47	OneQ-SL	Ø4.2X12mm	0	Cerasorb M Titanium Mesh

**Pre-operation**




**Fig1.** Pre-op radiography




**Fig2.** Pre-op intraoral photo


**Implant 1<sup>st</sup> surgery**



**Fig3.** After drilling, drilling path was checked by paralleling pines



**Fig4.** Drilling



**Fig5.** Implant Installation with OneQ-SL  $\varnothing$ 4.2X12mm



**Fig6.** Non-resorbable barrier membrane coverage



**Fig7.** Suture

**Implant 2<sup>nd</sup> surgery**



**Fig8.** Flap reflection for 2<sup>nd</sup> OP



**Fig9.** Healing abutment connection

**Prosthesis Delivery**



**Fig10.** Connected abutment



**Fig11.** Delivery prosthesis intraoral photo



**Fig12.** Delivery prosthesis radiography




Fig10. Delivery prosthesis CT

### Conclusion

Here, we attempted to find out if GBR with Mesh using horizontal augmentation was possible using only a synthetic bone graft. Most augmentations can be performed using any bone material with no concerns if placed under the mesh. However, many dentists still express concern using synthetic bone grafts. In this case, synthetic material was used, along with mesh and a DENTIS OneQ-SL Implant Maxillary posterior restoration was performed along with sinus elevation and simultaneous implant installation of the DENTIS OneQ-SL implant. No problems were noted following the final delivery of the prosthesis. This case exhibited successful sinus elevation with simultaneous implant installation. No concerns were noted after delivery of the prosthesis.

# Maxillary Posterior Restoration with Sinus Elevation and Simultaneous Implant Installation

Dae-Sung Kim  
 DDS, MSD.Ph.D, Director, MediPlant D.C / MEDI Academy of Implant Dentistry

## Patient Information

<b>Placement Implant Area</b>	6 7		<b>Age/Sex</b>	52Y/Male	
<b>C.C.</b>	Chronic periodontitis Loss of posterior teeth				
<b>Treatment Plan</b>	Sinus elevation and simultaneous implant installation				
<b>Material and Methods</b>	<b>Tooth No.</b>	<b>Product</b>	<b>Fixture Size</b>	<b>Sinus OP</b>	<b>Bone Graft Material</b>
	#26	OneQ-SL	Ø4.8X10mm	0	80% Cerasorb M +20% BioSS
	#27	OneQ-SL	Ø4.8X10mm	0	80% Cerasorb M +20% BioSS

**Pre-operation**




Fig1. Radiography before extraction



Fig2. Pre-op radiography

**Intra-operation**



**Fig3.** Sinus lateral window open for sinus augmentation



**Fig4.** Filled the bone material to the sinus after implant placement



**Fig5.** Cover with the cover screw



**Fig6.** Resorbable barrier membrane coverage



**Fig7.** Suture



**Post-operation**



Fig8. Post-op radiography



Fig9. Connected abutment



Fig10. Delivery prosthesis radiography

**Conclusion**

Maxillary posterior restoration by sinus elevation and simultaneous implant installation with DENTIS OneQ-SL. And still has no problem after delivery final prosthesis.

# Sinus Elevation and Simultaneous Implant Installation with OneQ-SL

Dae-Sung Kim  
 DDS, MSD.Ph.D, Director, MediPlant D.C/MEDI Academy of Implant Dentistry

## Patient Information

<b>Placement Implant Area</b>		6	<b>Age/Sex</b>	54Y/Male	
<b>C.C.</b>	Chronic periodontitis and bone destruction				
<b>Treatment Plan</b>	Sinus elevation and simultaneous implant installation				
<b>Material and Methods</b>	<b>Tooth No.</b>	<b>Product</b>	<b>Fixture Size</b>	<b>Sinus OP</b>	<b>Bone Graft Material</b>
	#26	OneQ-SL	Ø5.2X10mm	0	Cerasorb M

**Pre-operation**



Fig1. Radiography before extraction



Fig2. Pre-op radiography

**Intra-operation**



**Fig3.** Implant placement and outlining sinus lateral window



**Fig4.** Fill the bone material



**Fig5.** Resorbable barrier membrane coverage



**Fig6.** Suture

**Post-operation**



**Fig7.** Post-op radiography OneQ-SL Ø5.2X10mm

## Prosthesis Delivery



Fig8. Delivery prosthesis radiography



Fig9. Connected abutment



Fig10. Delivery prosthesis intraoral photo

## Conclusion

This is a sinus elevation and simultaneous implant installation case. After delivery prosthesis, there has no problem.

---

[www.dentimplant.co.kr](http://www.dentimplant.co.kr)

---



DENTIS Co., LTD.  
(Headquarters)

---

Tel. +82-53-583-2804~5  
Fax. +82-53-583-2806  
99, Seongseoseo-Ro, Dalseo-Gu, Daegu, Korea

---

DENTIS USA

---

Tel. +1-323-677-4363~5  
Fax. +1-323-677-4366  
6 Centerpointe Drive, Suite 600 La Palma CA 90623

---

DENTIS TAIWAN

---

Tel. +886-2-2808-5933  
Fax. +886-2-2624-2309  
2F.,No.48,31,Minzu Rd., Tamsui Dist, New Taipei City 251,Taiwan(R.O.C.)

---

DENTIS CHINA

---

Tel. +021-5111-3828  
Fax. +021-5111-3828  
上海市长宁区中山西路933号2205室

---