

For education

SIMPLE GUIDE **Plus**

User Manual Ver.1





CONTENTS

| | |
|---|----|
| I . Introduction | 4 |
| II . Advanced Preparation | 11 |
| Step1. CT Scan | |
| Step2. Scanning | |
| III . Method for chair-side fabrication | 14 |
| Step1. Implant Planning | |
| Step2. Surgical Stent 3D Printing | |
| IV . Surgical Procedure | 22 |
| Simple Guide Plus KIT | |
| Drilling Sequence | |
| Surgical Process | |
| V . Clinical Cases | 40 |
| Case1. Sleeveless Case | |
| Case2. Sinus OP Case | |
| Case3. Implant Placement Avoid Sinus Case | |
| Case4. Implant Placement Avoid Sinus Case | |
| Case5. Immediate Implant Placement Case | |
| Case6. Multiple Implant Placement Case | |
| VI . Specifications | 60 |



I . Introduction

What is **SIMPLE GUIDE Plus?**

Successful
Digital
Guide Surgery

Perform
a precise
surgery



✚ **SIMPLE GUIDE Plus** simple, convenient digital implant system with higher accuracy and compatible with all dental implant systems



Completion of 1 visit guide solution from diagnosis to guide output, and to surgery as well in dental clinic.

Economical price of the kit and cost-effective composition of the drills.

High accuracy without data loss due to digitalization of data of the entire processes of treatment protocol.

Guide solution in which clinical know-hows of digital clinicians and Dentis are concentrated.

Completion of open guide system due to compatibility of CEREC Guide 2

Easy guide planning system

SIMPLE GUIDE **Plus**

Why SIMPLE GUIDE Plus is Special?

Real One Day System

Only one-day guide system from diagnostic, surgical guide printing to implant surgery in your clinic.

- + Enable to surgery within 3hours(With model process) or 2hours(with intraoral scanner)
(*Single or couple implant case)
- + Depends on patient's condition, Immediate prosthetics loading (Coming soon)

Open Guide System

Open guide system compatible with all implant systems.

- + Compatible with CEREC Guide2
- + Make your own surgical guide in your clinic or guide center



Economic System

Accurate and cost-effective with convenient surgical kit.

- + Economic system with self guide production at chair side or center

Accurate System


The most accurate system with innovate design of the step drill & reamer drill.

- + 2 way drill and sleeve contacts improve accuracy
- + Controls space between guide and drill / sleeve and drill

User-friendly System

Convenient and easy system with new planning software.

- + Improve DACOM and STL merge algorithm to design surgical guide in 10 minutes
- + Convenient and easy software
- + Optional titanium sleeve and sleeve-less system
- + Optional open and close guide tube



Work-flow



CBCT Data



Scanning



Planning



3D Printing



Surgery



Temporary

Material



CT



Intraoral or Model scanner



PC & S/W



ZENITH U or Zenith D



Simple Guide PLUS
& Final Drill



ZENITH U or Zenith D

II . Advanced Preparation

Step1. CT Scan

Step2. Oral or Model Scanning

Step1. CT Scan

+ Material

| | |
|-------------------------------|--|
| Essential Material | CT, Cotton roll |
| Recommended Resolution | Slice interval < 0.25 mm Slice thickness < 0.25 mm Resolution < 0.25 mm x 0.25 mm |
| Recommended FOV Size | 10X8.5 or higher (When it is below that level, it is impossible to fabricate the guide up to full arch.) |



+ Method

1. CT Scan

- ① Remove any metallic objects in the head and neck region such as earrings, and glasses.
- ② Assure jaw separation during scan by biting over a cotton roll or similar object (Fig.1, 2, 3).



Fig.1 Open bite




Fig.2 Closed bite



Fig.3 Cotton roll

Notice If the patient moves during CT taking and the CT image is thus shaken, or CT FOV size is small (10 × 8.5 and over is recommended), or CT is cut out, or there are many prosthesis, it may be difficult to fabricate the guide.



Fig.4 Correct CT data




Fig.5 Shacked CT data




Fig.6 Unenough CT Information



Fig.7 Many Metallic Restorations Cases

2. Extract CT File (Dicom file)

Notice Method for file extraction may differ depending on the products. For detailed method, check on each manufacturers or inquire additionally.

Step2. Oral or Model Scanning

+ Material

| | |
|--------------------|--|
| Essential Material | (Option1) Intra-oral scanner |
| | (Option2) Stone model or model scanner |



+ Method

1. Intra-oral Scan Data



Extract STL File



2. Model Scan Data



Stone Model Scan



Extract STL File





III . Method for chair-side fabrication


Step1. Implant Planning

Step2. Surgical Stent 3D Printing

Step1. Implant Planning

+ Outline

- Guide planning S/W is a pre-operative planning software used to plan the placement of one or more implants based on CT and 3D scan data.
- The implant library that includes authenticated implants are supported.
- User can export the aligned implant and scan data as STL file.
- The use of this software requires having the necessary expertise in implant dentistry.



+ Recommended Guide Planning S/W

- DentiqGuide
- CEREC GALILEOS

Notice Method for file extraction may differ depending on the products. For detailed method, check on each manufacturers or inquire additionally.

Step1. Implant Planning

+ Material

| | |
|--------------------|---------|
| Essential Material | PC, S/W |
|--------------------|---------|

1. System Requirements

| Category | Minimum System Requirements | Recommended |
|---------------|------------------------------------|---|
| OS | Windows 7, 8, 10 (64bit) | Windows 7, 8, 10 (64bit) |
| RAM | 8GB | 16GB or higher |
| Graphics Card | cs card compatible with DirectX 11 | High end graphic card compatible with DirectX 11 or higher GPU-dedicated memory 2GB |
| CPU | Intel i3 Dual Core | Intel i5 Quad Core or higher |
| HDD | 2GB free space | 4GB free space |

* Minimum system requirements are only available when the S/W to be in sole occupancy.

+ Work-flow



Reference : DentiqGuide

Step2. Surgical Stent 3D Printing

+ Material

| | | |
|--------------------|-------------|--|
| Essential Material | Printing | ZENITH D or U, resin, vat sealer, vat, plate, PC, scraper, urethane paddle |
| | Cleaning | Cleaning container, alcohol, brush, ultrasonic waves |
| | Fabricating | UV curing, sleeve, handpiece |

1. ZENITH U

- Multi play desktop 3D printer that could be utilized in wide range of application



| | |
|------------------------|---|
| Technology | Stereolithography Apparatus (SLA) |
| Build Size | 110 x 110 x 150mm |
| Machine Size X/Y/Z | 354 x 366 x 483mm |
| Weight | 17.5 kg |
| Interface | USB |
| Layer Thickness | 16 μm, 50 μm, 100 μm |
| Light Source | Blue Laser |
| Software | ZENITH S/W |
| OS | Windows 7,8,10 / Mac OS 10.10 or higher |
| Input file format | Standard STL file |
| Power supply | Input : AC 100~240V / 50-60Hz, Output : DC 24V 5A |
| Material Properties | Photopolymer Resin |
| Electrical consumption | 120W |



2. ZENITH D

- Popular 3D printer special for model and guide
- DLP method 3D printer that is economical and cost effective



| | |
|---------------------|---|
| Technology | DLP(Digital Light Processing) |
| Build Size | 128 x 80 x 150 mm |
| Machine Size X/Y/Z | 340 x 460 x 430 mm |
| weight | 15 kg |
| XY Resolution | 100 μm |
| Interface | USB, Wi-Fi |
| Layer Thickness | 50 μm, 100 μm |
| Light Source | 405 μm LED |
| Software | Zenith D Slicer |
| OS | Windows 7,8,10 / Mac 10.10 or higher |
| Input file format | Standard STL file |
| Power supply | INPUT : AC 100V-240V / 50-60 Hz 2A max OUTPUT : 24V 5A |
| Material Properties | Photopolymer Resin |



Notice Specifications of the product may change depending on the circumstances of the company.

Step2. Surgical Stent 3D Printing

+ Method

1. 3D Printer S/W

- ① Import STL file by executing 3D printer S/W.
- ② Determining the location of the output, Designing the strut and floor surface strut.
- ③ Data transfer by using 3D printer(Automatic start-up)



2. Printing


① 3D Printer



② Removal of output board



③ Separation of output



3. Post-treatment

① Cleaning



② Separation of the strut



③ UV Curing



Notice For 3D printer manual, refer ZENITH software manual and user manual.

Step2. Surgical Stent 3D Printing

+ Guide Stent Fabrication

1. Sleeveless Guide Stent Trimming (Optional)

- It is used when the hole of the surgical stent output is not printed out precisely.
- It is used to lower the ratio of error while performing surgery by maintaining the gap between surgical stent and drill evenly.

① Prepare surgical stent, Ø4.5 Reamer of SIMPLE Guide Plus KIT



Surgical Stent



Ø 4.5 Reamer

② Perform drilling of the surgical stent under 1000–1200 rpm with Ø 4.5 Reamer extra-orally.

Notice * Perform reaming so that it becomes perpendicular to the guide.
* Pay attention to the path.



Step2. Surgical Stent 3D Printing

2. Sleeve Guide Stent Trimming (Optional)

- It is used when the hole of the surgical stent output is not printed out precisely.
- It is used to lower the ratio of error while performing surgery by maintaining the gap between surgical stent and drill evenly.

- ① Prepare Surgical guide stent, Ø5.3 Reamer of SIMPLE Guide Plus KIT의 Ø5.3 Reamer, Sleeve Connector, Sleeve



Surgical Guide Stent



Ø 5.3 Reamer



Sleeve (Inside Ø4.5/ Outside Ø5.3)



Sleeve Connector

- ② Perform drilling of the surgical stent under 1000–1200 rpm with Ø 5.3 reamer extra-orally.

Notice Perform reaming so that it becomes perpendicular to the guide. Pay attention to the path.



- ③ Connect the sleeve into the surgical stent by connecting the sleeve into the sleeve connector.

Notice Removing and attaching is easy without separate sleeve adhesive.



3. Sleeve Removing and Attaching Method

- ① Prepare sleeve connector
- ② Removal is easy by inserting the sleeve connector into the opposite side of the stent.





IV. Surgical Procedure

Simple Guide Plus KIT
Drilling Sequence
Surgical Process

1. SIMPLE GUIDE Plus KIT




Simple Guide Plus KIT

2. Component

· Reamer Drill (Option)

- Maintain the tolerance steadily between guide and drill, between sleeve and drill, and between sleeve and guide
- It is used when surgical stent is printed inaccuracy
- \varnothing 4.5 is used in sleeveless/ \varnothing 5.3 is used in sleeve
- Code : \varnothing 4.5_DSGPR45, \varnothing 5.3_DSGPR53




· Tissue Punch

- Remove soft tissue for flapless surgery
- Bone flattening function
- Code : \varnothing 4.0_DSGPTP




· Initial Drill

- Mark accurate points at any bone density
- Remove the gingiva that has not been removed by tissue punch
- Code : \varnothing 1.8_DSGPID



· Pilot Drill


- Enlarge depth and diameter
- Minimize bone cutting through drill sequence avoid bone heating
- Code : \varnothing 2.0x8.5 (16.5)_DSGPSD 2016
 \varnothing 2.0x10 (18)_D SGPSD 2018
 \varnothing 2.0x11.5 (19.5)_DSGPSD 2019



Simple Guide Plus KIT


• Straight Drill

- Step drill of which error ratio is dramatically reduced due to 2 point fixation
- Only with single final drilling, implant installation is possible
- Code : $\varnothing 2.8 \times 10.5$ (18.5)_DSGPSD 2818
 $\varnothing 2.8 \times 12$ (20)_DSGPSD 2820
 $\varnothing 2.8 \times 13.5$ (21.5)_DSGPSD 2821
 $\varnothing 3.5 \times 10.5$ (18.5)_DSGPSD 3518
 $\varnothing 3.5 \times 12$ (20)_DSGPSD 3520
 $\varnothing 3.5 \times 13.5$ (21.5)_DSGPSD 3521
 $\varnothing 4.3 \times 10.5$ (18.5)_DSGPSD 4318
 $\varnothing 4.3 \times 12$ (20)_DSGPSD 4320
 $\varnothing 4.3 \times 13.5$ (21.5)_DSGPSD 4321



• Fixture Driver (Machine / Ratchcet) (Optional)

- It is composed of 2 types which are for machine and for ratchet
- 3 points marking enable to implant hex position
- Enable to adjust depth with laser Marking
- When the fixture that is smaller than $\varnothing 4.2$ is installed, direct use on the stent is possible. If it is over that level, using is possible after removing the sleeve
- It is compatible with domestic implant fixtures (Choice is possible depending on using product)
- Code : Machine ($\varnothing 4.5$)_DSGPDA, Ratchet ($\varnothing 4.5$)_DSGPRDA
Machine ($\varnothing 4.5$)_DSGPDO, Ratchet ($\varnothing 4.5$)_DSGPRDO



• Profile Drill

- To trim residue bone after implant insertion for best supra structure parts connection
- Trim up to 6.5 (Regular implant)
- Code : $\varnothing 4.0$ _DSGPPD



Simple Guide Plus KIT

3. Separately sold product

- **Anchor Drill**

- A drill used for inserting the screw for fixation of surgical stent
- Ø1.3
- Code : DGATD13



- **Anchor Driver**

- A tool used for fixation of anchor screw
- Two different types for machine and for ratchet
- Code : Machine_DMHD24
Ratchet_DRHDS24



- **Anchor Screw**

- A tool used for fixation of surgical stent
- Code : DGAS18



- **Sleeve**

- Separate sleeve adhesive is unnecessary, and the groove on labial and lingual side makes checking hexa-surface easy
- Code : DSGPS



- **Sleeve Connector**

- A tool for removing and attaching the sleeve
- Removal is possible when insertion is done on the opposite side of the stent
- Code : DSGPSC



Simple Guide Plus KIT

4. Method for Storage and Management of the KIT

1) Cleaning and disinfection

- ① Thoroughly pre-rinse blood stain or foreign body on the instruments after using the kit by using a cleaning brush on the surface in distilled water or 30 ~ 40°C running water for 20 seconds.
- ② Pre-rinse it by immersing it in disinfectant liquid for 10 minutes.
- ③ Cleanse additionally by using ultrasonic cleanser.
- ④ After cleansing it by using detergent and cleaning brush, wash it in running water thoroughly.
- ⑤ Either 100% naturally dry the cleansed instruments or use a clean cloth to directly remove wetness.
- ⑥ Reposition the dried instruments in accordance with the mark of base plate of the kit.
- ⑦ Wrap the kit with sterilization wrap.
- ⑧ Mark the sterilization date after attaching sterilization tape on sterilization wrap.
- ⑨ Place the wrapped kit into the sterilization device and proceed sterilization.




2) Method for storage

- ① Store it in room temperature on uncontaminated area.
- ② Check the marked sterilization date, and if it has not been used within 3 ~ 4 days, re-sterilize it before using it for surgery.

Drilling Sequence

Universal System



Notice According to the patient's bone quality, over drilling or under drilling


Drilling Sequence

DENTIS OneQ Implant System


One System for All
Implant Indications

OneQ-SL

Ø3.0 x 8.0mm



Ø3.3 x 8.0mm




Drilling Sequence

DENTIS OneQ Implant System


One System for All
Implant Indications

OneQ-SL

Ø3.9 x 8.0mm



Ø4.2 x 8.0mm




Drilling Sequence

DENTIS OneQ Implant System


One System for All
Implant Indications

OneQ-SL

Ø4.7 x 8.0mm



Ø5.2 x 8.0mm



Surgical Procedure

1. Preparation of Guide Stent

1) Chlorohexidine disinfection of guide stent




- Prior to surgery, immerse the surgical stent in chlorohexidine for 5 – 10 minutes for disinfection.

Surgical Procedure

2. Remove Soft Tissue


1) Incision & Flap elevation (It could be chosen depending on the situation)



- When performing flap surgery, make an incision on the gingiva of the site for implant insertion by using a blade, and separate periosteum by using periosteal elevator or similar kind instrument.

Notice In flapless surgery case, use tissue punch.

2) Tissue Punch (It could be chosen depending on the situation)




- When performing flapless surgery, remove the gingiva on the site for implant installation by using tissue punch.
- Healing is faster since the wound is smaller than that by incision and flap elevation.

Tip Since tissues punch has bone flattening function as well, bone surface can be flattened without additional flattening drilling procedure

Notice Drilling is done under 800~1000 rpm recommendation

Surgical Procedure

3. Guide Stent Installation



Surgical Guide Stent

- Installation of the patient's surgical stent that has been fabricated.

Tip * Installation is done in planned direction on guide S/W image design. After installation, check whether there is opened part by pressing bilateral ends of the stent through open window, and nicely settle it in the mouth.

4. Drilling

Notice Drilling is done under 800~1000 rpm if it is high speed drilling, and 50~100 rpm if it is low speed drilling. While drilling is being done, insert the drill in the guide hole and check whether it has been accurately fixed.
All drillings must be done under copious water irrigation with pumping action (up and down).

1) Initial Drill Drilling




- It is a drill that is used for forming an accurate hole on the initial stage, and it could be used regardless of quality or form of the bone.

Notice

Drillings must be done under copious water irrigation with pumping action (up and down).

2) Pilot Drill Drilling




- Ø2.0 Pilot drill for widen diameter and depth

Notice

- * Appropriate size drill for fixture length is chosen, and drilling is not done up to the full size but to the size that is one step below to prevent bone heating.
- * Before drilling, adjust the drill into the pre-made hole and perform drilling with pumping action (up and down).

Surgical Procedure

5. Step Drill Drilling



Notice Before drilling, adjust the drill into the pre-made hole and perform drilling with pumping action (up and down).

6. Removal Guide Stent



· Removal Surgical stent.

Surgical Procedure


7. Final Drill Drilling



Final Drill

- Perform drilling in accordance with the method recommended by the manufacturer by using final drill of fixture surgical kit.

8. Countersink Drilling (Optional)




Countersink

- Perform drilling in accordance with the method recommended by the manufacturer by using countersink drill.

Surgical Procedure

9. Installation



- Fixture installation with fixture driver.
- When other implant system is used, however, and the fixture driver or separately purchased drivers are thus not compatible, use the driver of the corresponding manufacturer for installation.

Notice * When the implant with $\varnothing 4.2$ diameter or smaller is being installed, final installation will be possible when guide stent is connected again after final drilling.

* If you install over $\varnothing 4.2$ diameter fixture, remove the guide stent and install the implant using machined fixture driver first and try-in guide stent again and install with ratchet driver.


10. Connect upper prophetic part



- Choose between cover screw and healing abutment, and connect upper prophetic part.

Surgical Procedure

11. Profile Drilling (Optional)



- When the upper prophetic part cannot connect because of bone, cut the bone using profile drill.

Notice Drilling is done under 800~1000 rpm recommendation



V. Clinical Cases

Case1. Sleeveless Case

Case2. Sinus OP Case

Case3. Implant Placement Avoid Sinus Case

Case4. Implant Placement Avoid Sinus Case


Case5. Immediate Implant Placement Case

Case6. Multiple Implant Placement Case

Case1. Sleeveless Case

#37 Implant Placement

Dr. Sangjin Suh / Yein Dental Clinic

| Patient Information | | | |
|------------------------|---|-----------|--------------|
| Placement Implant Area |  | Sex / Age | Female / 39Y |
| C.C. | #37 tooth missing | | |
| Treatment Plan | #37 implant placement using surgical guide stent and GBR | | |
| Materials and Methods | <div>1. Guide design using DentiqGuide S/W</div> <div>2. Surgical guide stent printed by 3D printer, Zenith</div> <div>3. Try-in guide stent and drilling with Simple Guide Plus KIT</div> <div>4. Removed guide stent and drilling with final drill of OneQ surgical guide KIT</div> <div>5. OneQ Implant placement (4.7X8mm)</div> <div>6. 3 months after connect metal customized abutment and delivery prosthesis</div> | | |

Pre-operation



Fig.1 Pre-op panorama

Guide Stent Fabrication Procedure

DentiqGuide S/W Planning




Fig.2 Alignment




Fig.3 Crown placement




Fig.4 Implant placement




Fig.5 Guide

3D Printer Printing



Fig.6 Guide stent printed by ZENITH

Implant Installation



Fig.7 Drilling with initial drill



Fig.8 Drilling with pilot drill



Fig.9 Drilling with step drill



Fig.10 Drilling with final drill



Fig.11 4.7x8mm OneQ-SL implant installation

Post-operation



Fig.12 Post-op CT



Fig.13 Post-op panorama



Fig.14 3 Ms after surgery



Fig.15 Connect metal customized abutment



Fig.16 Delivery prosthesis

Conclusion

In the case of surgery on mandibular 2nd molar, there often are the cases where implant installation while securing the prosthetic position is difficult due to limitation of sight securing and mouth opening by anterior teeth. However, through simple guide plus system, guide surgery could be performed easily and precisely on the planned position.

Case2. Sinus OP Case

Upper Molar Implant Placement

Dr. Sangjin Suh / Yein Dental Clinic

| Patient Information | | | |
|------------------------|--|-----------|------------|
| Placement Implant Area | 7 6 4 | Sex / Age | Male / 55Y |
| C.C. | #16, #17 tooth missing, #14 root rest | | |
| Treatment Plan | #14, #16, #17 implant placement using surgical guide stent | | |
| Materials and Methods | 1. # 14 r.r extraction 2. Guide design using DentiqGuide S/W 3. Surgical guide stent printed by 3D printer, Zenith 4. Try-in guide stent and drilling with Simple Guide Plus KIT 5. Removed guide stent and drilling with final drill of OneQ surgical guide KIT 6. OneQ Implant placement #14, #17 (4.7X10mm), #16 (4.7X8mm) | | |

Pre-operation



Fig.1 Pre-op panorama



Fig.2 Pre-op clinical view

Guide Stent Fabrication Procedure

DentiqGuide S/W Planning



Fig.3a-c Implant placement




Fig.4a-b Guide

3D Printer Printing



Fig.5a-b Guide stent printed by ZENITH

Implant Installation



Fig.6a-b Drilling with initial drill



Fig.8 Drilling with step drill



Fig.9 Drilling with final drill

Fig.10 #14, #17 (4.7X10mm), #16 (4.7X8mm)
OneQ-SL implant installation and connected healing abutment and suturing

Post-operation



Fig.11 Post-op panorama



Fig.12 Post-op 6 months

Conclusion

It is a case of immediate implantation in sinus elevation by using simple guide plus system. Imaginary prosthesis was made and implant position was planned. By using metal sleeve, tolerance between drill and sleeve was reduced and thus stable drilling was possible. It was possible to install the implant on a precise location by using the guide after extracting the remaining root of right maxillary 1st premolar. And, on the area on right maxillary 1st pre-molar, implantation was planned up to lower border of maxillary sinus, and bone grafting was safely done after performing the drilling so that the mucosa of maxillary sinus is not ruptured by guide surgery. During sinus elevation procedure, the height of remaining bone that had been marked while implant planning could be referred.

Case3. Implant Placement Avoid Sinus Case

Upper Molar Implant Placement Avoid Sinus

Dr. Sangjin Suh | Yein Dental Clinic

| Patient Information | | | |
|------------------------|---|-----------|-----------|
| Placement Implant Area | 7 5 | Sex / Age | Male / 48 |
| C.C. | #15, #17 chronic periodontitis | | |
| Treatment Plan | #15, #17 implant placement using surgical guide stent | | |
| Materials and Methods | 1. # 15, #17 extraction 2. Guide design using DentiqGuide S/W 3. Surgical guide stent printed by 3D printer, Zenith 4. Try-in guide stent and drilling with Simple Guide Plus KIT 5. Removed guide stent and drilling with final drill of surgical guide KIT 6. Implant placement #15 (4.0X10mm), #17 (5.0X10mm) 7. Bone graft application 8. 3 months after connect metal customized abutment and delivery prosthesis | | |

Pre-operation



Fig.1 Before extraction panorama



Fig.2 Pre-op panorama



Fig.3 Pre-op CT



Fig.4 Pre-op clinical view

Guide Stent Fabrication Procedure

DentiqGuide S/W Planning



Fig.5a-b Implant placement



Fig.6 Guide

3D Printer Printing



Fig.7 Guide stent printed by ZENITH

Implant Installation



Fig.8 Try-in Guide Stent



Fig.9 Drilling and implant placement



Fig.10 Connect healing abutment and GBR



Fig.11 Suturing

Post-operation



Fig.12a-b Post-op CT. a. #16, b. #17



Fig.13 3 Ms after surgery



Fig.14 Connect metal customized abutment



Fig.15 Delivery prosthesis (Full zirconia)

Conclusion

In this case, although implantation was planned on #15 to 17 3-unit bridge after extraction due to chronic periodontitis, elevation of inferior border of maxillary sinus and ridge augmentation had to be performed and long term healing period over 6 months had to be waited.

In this case, however, tilted implant was planned by calculating the possible angle that allows custom made abutment on posterior portion where remaining bone is preserved by using surgery guide, and quantity of bone graft was thus minimized, and only simple ridge augmentation was done and the prosthesis was completed in 3 months.

Case4. Implant Placement Avoid Sinus Case

Upper Molar Implant Placement Avoid Sinus

Dr. Sangjin Suh | Yein Dental Clinic

| Patient Information | | | |
|------------------------|---|-----------|------------|
| Placement Implant Area | <div><div></div><div>67</div></div> | Sex / Age | Male / 69Y |
| C.C. | #26, #27 implant placement using surgical guide stent | | |
| Materials and Methods | <div><div>1. Guide design using DentiqGuide S/W</div><div>2. Surgical guide stent printed by 3D printer, Zenith</div><div>3. Try-in guide stent and drilling with Simple Guide Plus KIT</div><div>4. Removed guide stent and drilling with final drill of surgical guide KIT</div><div>5. Implant placement OneQ-SL #26 (4.7X10mm), #27 (5.2X8mm)</div></div> | | |

Pre-operation



Fig.1 Pre-op panorama



Fig.2 Pre-op clinical view

Guide Stent Fabrication Procedure

DentiqGuide S/W Planning




Fig.3 Alignment




Fig.4 Crown placement




Fig.5a-b Implant placement





Fig.6 Guide

3D Printer Printing



Fig.7a-b Guide stent printed by ZENITH

Implant Installation



Fig.8 Try-in Guide Stent



Fig.9 Implant placement

Post-operation




Fig.10a-b Post-op CT. a. #16, b. #17



Fig.11 Post-op panorama



Fig.12 Post-op 3 months panorama



Fig.13 3 months after metal custom abutment



Fig.14 Delivery final prosthesis


Conclusion

On #26 and #27, vertical height of the remaining bone is insufficient. In this case, conventional sinus elevation technique should be used. However, 10 mm and 8 mm length implants were able to be installed without bone grafting by positioning the implant being tilted toward the portion where bone quantity was sufficient by using the guide, and the surgery was thus done simply.

Case5. Immediate Implant Placement Case

#36 Extraction and immediate implant placement

Dr. Sangjin Suh | Yein Dental Clinic

| Patient Information | | | |
|------------------------|---|-----------|--------------|
| Placement Implant Area |  | Sex / Age | Female / 52Y |
| C.C. | #36 tooth fracture | | |
| Treatment Plan | #36 extraction, immediate implant placement using surgical guide stent and GBR | | |
| Materials and Methods | <ol style="list-style-type: none">1. Guide design using DentiqGuide S/W2. Surgical guide stent printed by 3D printer, Zenith3. Try-in guide stent and drilling with Simple Guide Plus KIT4. Removed guide stent and drilling with final drill of surgical guide KIT5. Implant placement (5.0X10)6. Bone graft and colla tape application | | |

Pre-operation




Fig.1 Pre-op panorama




Fig.2a-b Pre-op CT



Fig.3 Pre-op clinical view

Guide Stent Fabrication Procedure

DentiqGuide S/W Planning




Fig.4 Alignment




Fig.5 Tooth removal which will be tooth extraction




Fig.6 Crown placement



Fig.7 Draw curve




Fig.8 Nerve setting



Fig.9 Implant placement




Fig.10 Sleeve select




Fig.11 Guide Window open

3D Printer Printing



Fig.12a-b Guide stent printed by ZENITH

Implant Installation



Fig.13 #36 extraction



Fig.14 Try-in Guide Stent



Fig.15 Drilling with initial drill



Fig.16 Drilling with pilot drill



Fig.17 Drilling with step drill



Fig.18 Drilling with final drill



Fig.19 Implant installation



Fig.20 GBR



Fig.21 Collar tape

Post-operation



Fig.22 Post-op CT



Fig.23 Post-op panorama



Fig.24 2 Ws after surgery



Fig.25 Post-op 2 months panorama



Fig.26 Delivery final prosthesis

Conclusion

It is an immediate implantation case of right mandibular 1st molar by using simple guide plus system. In immediate implant cases of mandibular molars, gaining stability of implant while drilling the septal bone accurately is often difficult. However, guide surgery could solve this issue. After installing implant, extraction socket was grafted with heterogenous bone, and soft tissue defect was closed with absorbable collagen membrane, and closure was done by using tissue adhesive. In 2 weeks, stable healing appearance is observed. The best part of simple guide plus is that the clinician himself / herself could directly plan surgery procedures and the outcome of the surgery will be satisfactory as well accordingly.

Case6. Multiple Implant Placement Case

Multiple implant placement

Dr. Sangjin Suh | Yein Dental Clinic

| Patient Information | | | | |
|------------------------|--|-----------|------------|--|
| Placement Implant Area | 6 4 2 2 6 6 4 5 6 | Sex / Age | Male / 70Y | |
| C.C. | #12, 14, 16, 22, 26, 34, 35, 36, 46 missing tooth | | | |
| Treatment Plan | Implant placement using surgical guide stent | | | |
| Materials and Methods | 1. Guide design using DentiqGuide S/W 2. Surgical guide stent printed by 3D printer, Zenith 3. Try-in guide stent and drilling with Simple Guide Plus KIT 4. Removed guide stent and drilling with final drill of surgical guide KIT 5. Implant Fixture #12, 22 (3.6X10), #14, 34, 35 (4.0X10), #16, 26 (5.0X10), #36, 46 (4.5X10) 6. Delivery prosthesis | | | |

Pre-operation



Fig.1 Pre-op panorama




Fig.2 Implant placement treatment plan

Guide Stent Fabrication Procedure

DentiqGuide S/W Planning




Fig.3a-d Upper Implant placement





Fig.4a-d Lower Implant placement




Fig.5 Upper Guide



Fig.6 Lower Guide

3D Printer Printing



Fig.7a-b Upper Guide stent printed by ZENITH



Fig.8a-b Lower Guide stent printed by ZENITH



Implant Installation



Fig.9 Try-in upper Guide Stent



Fig.10 Drilling



Fig.11 Upper Implant installation



Fig.12 Try-in lower Guide Stent



Fig.13 Drilling



Fig.14 Lower Implant installation

Post-operation



Fig.15 Post-op CT



Fig.16 Post-op panorama



Fig.17 Post-op clinical view

Conclusion


This is a case where it is not easy to install implant on a precise spot in prosthetic perspective without using the guide because number of remaining teeth is insufficient. Since #13 is mesially tilted and it needs to be extracted due to chronic periodontitis, imaginary extraction was done on guide software and then guide for surgery was fabricated. In order to check the condition of irregular alveolar bone, flap opening was done and then implant was installed by using surgery guide. Upper and lower implants were simultaneously installed, and the surgery could be finished quickly on the precise location thanks to the guide.



VI. Specifications

* SIMPLE GUIDE Plus Kits

- Exclusive Surgical Guide Kits of Dentis Simple Guide Plus System



| Type | Size | Code |
|------------|------------------------|-------|
| Guide Kits | 134(W X39(H)X109(D) mm | DSGPK |

* SG SCAN M

- Stable Scan Stage (SSS) method with moving camera technology allows you scan without any fixture or adhesive.



| Type | Size | Code |
|---------------|-------------------------|------------|
| Model Scanner | 330(W)X495(H)X430(D) mm | MEQDOF0001 |

* ZENITH U

- Independent development by DENTIS
- 16 μ m resolution galvanometer scanning method of 3D printer for specialized for the dental industry
- Grade optical system application
- Working Area: 110X110X150 (X,Y,Z/mm)



| Type | Size | Code |
|------|---------------|--------------|
| SLA | 354X366X483mm | PP1000-02120 |

* ZENITH D

- Independent development by DENTIS
- Self-developed ZENITH resins enable more reliable, low-cost production of a variety of dental products in-clinic, including surgical guides, dental models, temporary C&B.
- Working Area: 128 x 80 x 145mm (X,Y,Z/mm)



| Type | Size | Code |
|------|---------------|--------------|
| DLP | 340x460x430mm | PP2000-02000 |

*** ZENITH Accessory****· ZENITH U VAT**

- Application of special resin tank plate make possible to use as a semi-permanently.



| Type | Code |
|-----------|--------------|
| Vat tank | PP1000-0V120 |
| Vat Plate | PP1000-0V110 |

· ZENITH D VAT

- More than 20 ~ 30 times use for printing



| Type | Code |
|----------|--------------|
| Vat tank | SP2MVA-10010 |

*** ZENITH Resin****ZMD-1000B Clear SG****· For ZENITH U**

- Working time (Min~Max) : 150 min. ~ 200 min. (100 μ m, Standard 2 full arch)
- The number of output arrangement on the working plate (Max) : 2-3 Full arch
- Each consumption (g) : About 15g (standard 1 full arch)

· For ZENITH D

- Working time (Min~Max) : 40 min. ~ 80 min. (100 μ m, Standard 2 full arch)
- The number of output arrangement on the working plate (Max) : 2-3 Full arch
- Each consumption (g) : About 15g (standard 1 full arch)



| Application | Weight | Code |
|----------------|--------|--------------|
| Surgical Guide | 1Kg | PP0R1Z-45122 |

ZMD-1000B Temporary**· For ZENITH U**

- Working time (Min~Max) : 150 min. ~ 200 min. (100 μ m, Standard 25 single crown)
- The number of output arrangement on the working plate (Max) : 20-25 single crown (Standard temporary exclusive plate)
- Each consumption (g) : About 1-1.5g (standard 3 unit bridge, include support)

· For ZENITH D

- Working time (Min~Max) : 25 min. ~ 60 min. (Influenced by crown length)
- The number of output arrangement on the working plate (Max) : 20-25 single crown (Standard temporary exclusive plate)
- Each consumption (g) : About 1-1.5g (standard 3 unit bridge, include support)



| Application | Weight | Code |
|-------------|--------|--------------|
| Temporary | 1Kg | PP0R3Z-45141 |

www.dentisimplant.co.kr



Tel. +82-1899-2804 | Fax. +82-53-583-2806

99, Seongseoseo-Ro, Dalseo-Gu, Daegu, Korea



Tel. +1-323-677-4363 | Fax. +1-323-677-4366

6 Centerpointe Drive, Suite 600 La Palma CA 90623



Tel. +021-5111-3828 | Fax. +021-5111-3828

上海市长宁区中山西路933号2205室

Specifications are subject to change without notice.

Copyright © DENTIS Co., LTD.

D-M-SGP-V1-201807-ENG